

logo not found or type unknown

Firma: Player Plus doo
Adresa: Svetogorska 9
Telefon: +381 11 3347 442
Fax: +381 11 3347 615
PIB: 106966344
E-mail: porudzbine@player.rs

Goldring E3 Cartridge Gramofonska Glava

Goldring E3 Cartridge Gramofonska Glava

Šifra: 14482
Kategorija proizvoda: Igle - Glave i Oprema
Proizvođač: Goldring

Cena: **17.880,00** rsd

Magnetic Duplex Technology™

All cartridges in the range feature Magnetic Duplex Technology™ -a 'dual magnet' generator assembly with an identical geometry to that of the original cutting head. This means it can more accurately reproduce the stereo information in the record groove, providing better channel separation and a more involving and musical sound stage.

Expert Adjustment

The E-series cartridges have been designed to allow micro-adjustments to be made during production using a special adjustment screw. This means the cartridge can be safely installed in a typical turntable because the angle that the stylus makes with the record surface has been pre-set. The result is an optimised channel balance with zero degrees of horizontal tracking adjustment required. This ensures a consistent and reassuring performance for the life of the product.

Stylus/Cantilever

Featuring spherical styli with a choice of cantilevers, the E1 and E2 models are ideal as an upgrade to any budget, pre-fitted or obsolete cartridge. The E3 model with its re-fined super elliptical stylus is a worthwhile upgrade for those wishing to squeeze even more information and detail from the record groove.



Colour

Violet

Cantilever

Aluminium

Stylus profile

Elliptical

Stylus radius

0.3 x 0.7 mil

Frequency range

20 Hz - 20 kHz

Channel balance

1.5 dB at 1 kHz

Channel separation

>20 dB at 1 kHz

Sensitivity

3.5 mV (1 kHz @ 5 cm/sec)

Static compliance

20 mm/N

Vertical tracking angle

20°

Load resistance

47 kΩ

Load capacitance

100 - 200 pF

Internal inductance

400 mH

Internal resistance

410 Ω

Cartridge mass

6.9 g

Fixing centres

12.7 mm (0.5 in)

Playing weight

1.5 g - 2.5 g (2.0 g recommended)